

BACTERICIDAL

In the presence of 5% blood serum and a contact time of 2 minutes, is effective against:

Pseudomonas aeruginosa

Salmonella enterica

Staphylococcus aureus

In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:

Carbapenem-resistant Klebsiella pneumoniae

Enterococcus faecalis Vancomycin Resistant [VRE]

Staphylococcus aureus - MRSA

Escherichia coli O157:H7

Listeria monocytogenes

Staphylococcus epidermidis - MRSE

In the presence of 5% blood serum and a contact time of 10 minutes, is effective against:

Acinetobacter baumannii

TUBERCULOCIDAL

Effective against *Mycobacterium bovis* (BCG) in 5 minutes at 20°C / 68°F on previously cleaned hard non-porous inanimate surfaces.

FUNGICIDAL

In the presence of 5% blood serum and a contact time of 2 minutes, is effective against:

Trichophyton interdigitale

In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:

Candida albicans

VIRUCIDAL*

In the presence of 5% blood serum and a contact time of 30 seconds, is effective against:

SARS-CoV-2 cause of COVID-19**

In the presence of 5% blood serum and a contact time of 3 minutes, is effective against:

Rhinovirus

Hepatitis B Virus

Hepatitis C Virus

In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:

Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)

Influenza A Virus [flu virus]

Respiratory Syncytial Virus [RSV]

Rotavirus

Feline Calicivirus [surrogate for Norovirus]

Adenovirus Type 2

Herpes Simplex Virus Type 1

Herpes Simplex Virus Type 2

Vaccinia Virus

In the presence of 5% blood serum and a contact time of 10 minutes, is effective against:

Canine Parvovirus

Poliovirus Type 1

DISPENSER DIRECTIONS

Open lid. Pull out wipe from center of roll and through lid. Close lid to retain moisture.

Item No. 490908



7 8 9 7 0 3 0 0 0 1 1

For non-emergency (e.g., current product information), call: 1-800-621-8829

Manufactured for: Preserve International

944 Nandino Blvd. Lexington KY 40511 U.S.A.

Preserve International is a wholly-owned subsidiary of Neogen Corporation
L4839-0922

COMPANION™ DISINFECTANT WIPES

- Plant-derived citric acid disinfectant
- Disinfects in as little as 30 seconds**
- Disinfects & deodorizes²

²when used according to disinfection directions

Active Ingredient:

Citric Acid.....0.60%

Other Ingredients.....99.40%

Total.....100.00%

EPA Reg. No. 34810-36-66171

EPA Est. No. 10088-WI-2

EPA Est. No. 92811-TN-1

EPA Est. No. 66171-TN-001

Beginning of Lot indicates EPA Est. No.

KEEP OUT OF REACH OF CHILDREN



DO NOT FLUSH



160 Pre-moistened Towelettes 6" x 7"

NET CONTENTS: 1 lb., 3 oz.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

TO CLEAN/DEODORIZE

Unfold premoistened cloth and wipe area to be cleaned for 30 seconds or until clean. Spot check to test for compatibility.

TO DISINFECT

Product is ready to use. Preclean areas with visible soil prior to disinfection. Unfold premoistened cloth and wipe hard non-porous area or item to be disinfected. Wet all surfaces thoroughly. Leave wet for 10 minutes. Wipe or let air dry. No rinsing or wiping is required except on food contact surfaces which require a potable water rinse after disinfection.

TYPES OF SURFACES/USE SITES

Suitable for use on washable, hard, non-porous surfaces of:

Storage Areas	Elevators
Home	Equine Farms
Cars	Exam Rooms
Airports	Grooming Facilities
Animal Care Facilities	Kennels
Animal Hospitals	Lobbies
Animal Labs	Office Breakrooms
Animal Life Science Laboratories	Offices
Bathrooms	Pet Animal Quarters
Break Rooms	Pet Areas
Breeding & Grooming Establishments	Pet Shops
Butcher Shops	Reception Areas
Cafeterias	Shelters
Central Supply Rooms	Trucks
Clinics	Vehicle Interiors
Delivery Trucks	Veterinary Clinics
Dental Offices	

Suitable for use on hard, non-porous surfaces:

Animal Equipment	Countertops
Automatic Feeders	Doorknobs
Cages	Exam Tables
Chairs	External Surfaces of medical equipment
Computer Surfaces	Litter Boxes

Medical Equipment surfaces/devices
Microscopes
MRI Equipment Surfaces
Pet Habitats
Pet Soiled Surfaces
Pet Toys

Stainless Steel
Stalls
Surgical Equipment Surfaces
Walls
Work Stations
X-ray equipment/surfaces

DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION: This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that [1] is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semicritical medical devices prior to sterilization or high-level disinfection.

***COMPANION DISINFECTANT WIPES KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS virus).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). PERSONAL PROTECTION: When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: Do Not Dilute. Leave precleaned surfaces wet for five minutes using mop, cloth or sponge to mitigate HIV-1 virus. Use a 5 to 10 minute contact time to kill all bacteria and viruses* in accordance with the organisms listed on the label.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

STORAGE AND DISPOSAL: Do not leave containers open during storage. No special handling is required. **Towelette:** Do not reuse towelette. Dispose of used towelette in trash. **Dispenser:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available.